

SCHEMA TECNICA

REAGENTI

Descrizione Reagenti		Quantità
1	Derivatizzante BSTFA + TMCS, 99:1 *20X1 ml	7
2	Std interno : Tricarballic acid 25 gr.	1
3	Std Interno: 2-Phenylbutyric acid 100 gr.	1
4	Hydroxylamine hydrochloride 500 gr. Reagent Plus 99% CAS n. 5470-11-1	1
5	Sodium Chloride 1kg ACS reagent 99% CAS n.7647-14-5	1
6	Piridina > 99,5% 1000 ml	1
7	Acido cloridrico 35% 1 Lt	2